## Enterolert\* Test Kit

**Introduction**


Enterolert\* detects enterococci, such as *E. faecium* and *E. faecalis*, in fresh and marine water. It is based on IDEXX’s patented Defined Substrate Technology\* (DST\*). When enterococci utilize their β-glucosidase enzyme to metabolize Enterolert’s nutrient-indicator, 4-methyl-umbelliferyl β-D-glucoside, the sample fluoresces. Enterolert detects enterococci at 1 cfu per 100 mL sample within 24 hours.

**Storage**

Store at 2–30°C away from light.

**Presence/Absence (P/A) Procedure**

- Add contents of one pack to a 100 mL sample in a sterile, transparent, nonfluorescing vessel.
- Cap vessel and shake.
- Incubate at 41 ±0.5°C for 24 hours.
- Read results according to Result Interpretation table below.

**Quanti-Tray\* Enumeration Procedure**

- Add contents of one pack to a 100 mL water sample in a sterile vessel.
- Cap vessel and shake until dissolved.
- Pour sample/reagent mixture into a Quanti-Tray\* or Quanti-Tray\*/2000 and seal in an IDEXX Quanti-Tray\* Sealer.
- Place the sealed tray in a 41 ±0.5°C incubator for 24 hours.
- Read results according to the Result Interpretation table below. Count the number of positive wells and refer to the MPN table provided with the trays to obtain a Most Probable Number.

**Result Interpretation**

Appearance	Result
Lack of fluorescence	Negative for enterococci
Blue fluorescence	Positive for enterococci

- Look for fluorescence with a 6–watt, 365 nm, UV light within 5 inches of the sample in a dark environment. Face light away from your eyes and towards the sample.
- Enterolert results are definitive at 24–28 hours. In addition, positives for enterococci observed before 24 hours and negatives observed after 28 hours are also valid.

**Procedural Notes**

- Marine water samples must be diluted at least tenfold with sterile fresh water. Multiply the MPN value by the dilution factor to obtain the proper quantitative result.
- Use only sterile, nonbuffered, oxidant-free water for dilutions.
- For comparison, a water blank can be used when interpreting results.
- This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures.
- Enterolert can be run in any multiple tube format. *Standard Methods for the Examination of Water and Wastewater*<sup>1</sup> MPN tables should be used to find Most Probable Numbers.
- Enterolert is a primary water test. Enterolert performance characteristics do not apply to samples altered by any pre-enrichment or concentration.
- Aseptic technique should always be followed when using Enterolert. Dispose of in accordance with Good Laboratory Practices.

**Quality Control Procedures**

- One of the following quality control procedures is recommended for each lot of Enterolert:
  - IDEXX-QC Enterococci<sup>2</sup>: *Enterococcus faecalis*, *Escherichia coli*, and *Streptococcus bovis*.
  - For each of the American Type Culture Collection (ATCC)<sup>3</sup> bacterial strains (*Enterococcus faecium* ATCC 35667, *Serratia marcescens* ATCC 43862, and *Aerococcus veridans* ATCC 10400), streak the culture onto labeled TSA or Blood Agar plates and incubate at 35±2°C for 18–24 hours.
    - For each bacterial strain, touch a sterile 1 µL inoculating loop to a colony and use it to inoculate a labeled test tube containing 5 mL of sterile deionized water. Close cap and shake thoroughly.
    - For each bacterial strain, take a 1 µL loop from the test tube and use it to inoculate a labeled vessel containing 100 mL of sterile deionized water. These are your controls.
- Follow the P/A Procedure or Quanti-Tray Enumeration Procedure above.
- Results should match the Result Interpretation table above.

**NOTE:** IDEXX internal quality control testing is performed in accordance with ISO 11133:2014. Quality Control Certificates are available at [idexx.com/water](http://idexx.com/water).

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1. Eaton AD, Clesceri LS, Rice EW, Greenberg AE, Franson MAH, eds. Standard Methods for the Examination of Water and Wastewater. 21st ed. (Standardverfahren für die Wasser- und Abwas-seruntersuchung). Washington, DC: American Public Health Association; 2005.

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Patentinformation: idexx.com/patents.

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<sup>[</sup>\*Enterolert, Defined Substrate Technology, DST und Quanti-Tray are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

<sup>[</sup>Patent information: idexx.com/patents.

<sup>[</sup>© 2015 IDEXX Laboratories, Inc. All rights reserved.

IDEXX Water Quality Control Laboratory is accredited to ISO/IEC 17025:2005

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IDEXX Laboratories, Inc., One IDEXX Drive, Westbrook, Maine 04092 USA